



DNASE2B Antibody, Biotin conjugated

Product Code	CSB-PA841227LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WZ79
Immunogen	Recombinant Human Deoxyribonuclease-2-beta protein (68-281AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Deoxyribonuclease-2-beta (EC 3.1.22.1) (DNase II-like acid DNase) (DNase2-like acid DNase) (Deoxyribonuclease II beta) (DNase II beta) (Endonuclease DLAD), DNASE2B, DLAD
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	DNASE2B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.