



DNASE2B Antibody, Biotin conjugated

Product Code	CSB-PA841227LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WZ79
Immunogen	Recombinant Human Deoxyribonuclease-2-beta protein (68-281AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Deoxyribonuclease-2-beta (EC 3.1.22.1) (DNase II-like acid DNase) (DNase2-like acid DNase) (Deoxyribonuclease II beta) (DNase II beta) (Endonuclease DLAD), DNASE2B, DLAD
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	DNASE2B