




Recombinant Human Interleukin-15 (IL15) (Active) (GMP)

Product Code	CSB-AP005861HU
Uniprot No.	P40933
Form	Lyophilized powder
Storage Buffer	Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4.
Product Type	Interleukin
Immunogen Species	Homo sapiens (Human)
Biological Activity	Fully biologically active when compared to standard. The ED50 as determined by a cell proliferation assay using murine CTLL-2 cells is less than 0.5 ng/ml, corresponding to a specific activity of > 2.0 × 10 ⁶ IU/mg.
Purity	> 98 % by SDS-PAGE and HPLC analyses.
Sequence	NWVNVISDLKKIEDLIQSMHIDATLYTESDVHPSCKVTAMKCFLLLELQVISLESG DASIHDTVENLILANNLSLSSNGNVTESGCKECELEEKNIKEFLQSFVHIVQMFI NTS
Research Area	Immunology
Source	E.Coli
Target Names	IL15
Expression Region	49-162aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag-Free
Mol. Weight	12.9 kDa
Protein Length	Full Length of Mature Protein
PubMed ID	8178155; 8932977; 8668345; 9405632; 14702039; 15489334; 21167836

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Endotoxin	Less than 0.01 EU/µg of rHuIL-15 GMP as determined by LAL method.
------------------	---

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.