



CCL21 Antibody, HRP conjugated

Product Code	CSB-PA004785DB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O00585
Immunogen	Recombinant Human C-C motif chemokine 21 protein (24-134AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	C-C motif chemokine 21 (6Ckine) (Beta-chemokine exodus-2) (Secondary lymphoid-tissue chemokine) (SLC) (Small-inducible cytokine A21), CCL21, SCYA21
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CCL21
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.