



TNNI3 Antibody, Biotin conjugated

Product Code	CSB-PA07939D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P19429
Immunogen	Recombinant Human Troponin I, cardiac muscle protein (2-210AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Troponin I, cardiac muscle (Cardiac troponin I), TNNI3, TNNC1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TNNI3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.