



RUNX1 Antibody, HRP conjugated

Product Code	CSB-PA020592LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q01196
Immunogen	Recombinant Human Runt-related transcription factor 1 protein (211-445AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Runt-related transcription factor 1 (Acute myeloid leukemia 1 protein) (Core-binding factor subunit alpha-2) (CBF-alpha-2) (Oncogene AML-1) (Polyomavirus enhancer-binding protein 2 alpha B subunit) (PEA2-alpha B) (PEBP2-alpha B) (SL3-3 enhancer factor 1 alpha B subunit) (SL3/AKV core-binding factor alpha B subunit), RUNX1, AML1 CBFA2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	RUNX1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.