




Recombinant Human Polyadenylate-binding protein 4 (PABPC4)

Product Code	CSB-EP615676HUe1
Abbreviation	Recombinant Human PABPC4 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q13310
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	MNAAASSYPMASLYVGDHSDVTEAMLYEKFSPAGPVLSIRVCRDMITRRSLG YAYVNFQQPADAERALDTMNFVVIKGGKPIRIMWSQRDPSLRKSGVGNVFIKNL DKSIDNKALYDTFSAFGNILSCKVVCDENGSKGYAFVHFETQEAADKAIEKMN GMLLNDRKVFVGRFKSRKEREAEELGAKAKEFTNVYIKNFGEEDDESLEKELFS QFGKTL SVKVMRDPNGKSKGFGFVSYEKHEDANKAVEEMNGKEISGKIIFVGR AQKKVERQAE LKRKFEQLKQERISRYQGVNLYIKNLDDTIDDEKLRKEFSPFGS ITSAKVMLEDGRSKGFGFVCFSSPEEATKAVTEMNGRIVGSKPLYVALAQRKE ERKAHLTNQYMQRVAGMRALPANAILNQFQPAAGGYFVPAVPQAQGRPPYY TPNQLAQMRPNPRWQQGGRPQGFQGMPSAIRQSGPRPTLRHLAPTGSECP DRLAMDFGGAGAAQQGLTDSCQSGGVPTAVQNLAPRAAVAAAAPRAVAPYK YASSVRSPHPAIQPLQAPQPAVHVQGGQEPLTASMLAAAPPQEQQKQMLGERLF PLIQTMHSNLAGKITGMILLEIDNSELLHMLESPESLRSKVDEAVAVLQAHHAKK EAAQKVGAVAAATS
Research Area	Others
Source	E.coli
Target Names	PABPC4
Expression Region	1-644aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag-Free
Mol. Weight	70.8 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.