




# Recombinant Human Syntaxin-binding protein 1 (STXBP1)

<b>Product Code</b>	CSB-EP022905HU
<b>Abbreviation</b>	Recombinant Human STXBP1 protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P61764
<b>Storage Buffer</b>	Tris-based buffer,50% glycerol
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥ 90% as determined by SDS-PAGE.
<b>Sequence</b>	MAPIGLKAVVGEKIMHDVIKKVKKKGEWKVLVVDQLSMRMLSSCCKMTDIMTE GITIVEDINKRREPLPSLEAVYLITPSEKSVHSLISDFKDPPTAKYRAAHVFFTD CPDALFNELVKSRAAKVIKTLTEINIAFLPYESQVYSLDSADSFQSFYSPHKAQM KNPILERLAEQIATLCATLKEYPAVRYRGEYKDNALLAQLIQDKLDAYKADDPT MGEGPDKARSQLLILDRGFDPSVPLHELTFQAMSYDLLPIENDVYKYETSGIG EARVKEVLLDEDDDLWIALRHKHIAEVSQEVTRSLKDFSSSKRMNTGEKTTMR DLSQMLKKMPQYQKELSKYSTHLHLAEDCMKHYQGTVDKLCRVEQDLAMGT DAEGEKIKDPMRAIVPILLDANVSTYDKIRIILLYIFLKNGITEENLNKLIQHAQIPP EDSEIITNMAHLGVPIVTDSTLRRRSKPERKERISEQTYQLSRWTPPIKDIMEDI EDKLDTKHYPYISTRSSASFSTTAVSARYGHWKHNKAPGEYRSGPRLIIFILGG VSLNEMRCAYEVTQANGKWEVLIGSTHILTPQKLLDTLKLNKTDEEISS
<b>Research Area</b>	Cancer
<b>Source</b>	E.coli
<b>Target Names</b>	STXBP1
<b>Protein Names</b>	Recommended name: Syntaxin-binding protein 1 Alternative name(s): MUNC18-1 N-Sec1 Protein unc-18 homolog 1 Short name= Unc18-1 Protein unc-18 homolog A Short name= Unc-18A p67
<b>Expression Region</b>	1-594aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-MBP-tagged
<b>Mol. Weight</b>	111.3 kDa
<b>Protein Length</b>	Full Length
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.