




Recombinant Chicken Antimicrobial peptide CHP1

Product Code	CSB-EP301330CH
Abbreviation	Recombinant Chicken Antimicrobial peptide CHP1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P80389
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Gallus gallus (Chicken)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	GRKSDCFRKSGFCAFLKCPSLTLISGKCSR FYLCCKRIR
Research Area	Others
Source	E.coli
Target Names	N/A
Protein Names	Recommended name: Antimicrobial peptide CHP1 Alternative name(s): Chicken heterophil peptide 1
Expression Region	1-39aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-KSI-tagged
Mol. Weight	19.8 kDa
Protein Length	Full Length

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.