




Recombinant Human PITH domain-containing protein 1 (PITHD1)

Product Code	CSB-EP867109HU
Abbreviation	Recombinant Human PITHD1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9GZP4
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 90% as determined by SDS-PAGE.
Sequence	MSHGHSHGGGGCRCAAEREPEPQGRGLAYGLYLRLIDLERLQCLNESREGSG RGVFKPWEERTDRSKFVESDADEELLFNIPFTGNVCLKGIIIMGEDDDSHPSSEM RLYKNIPQMSFDDTEREPDQTFSLNRDLTGELEYATKISRFSNVYHLSIHISKNF GADTTKVFYIGLRGEWTELRRHEVTICNYEASANPADHRVHQVTPQTHFIS
Research Area	Others
Source	E.coli
Target Names	PITHD1
Protein Names	Recommended name: PITH domain-containing protein 1
Expression Region	1-211aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	31.6 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.