




Recombinant Shigella flexneri Protein mxiH (mxiH)

Product Code	CSB-EP363183SZB
Abbreviation	Recombinant Shigella flexneri mxiH protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0A223
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Shigella flexneri
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MSVTVPNDDWTLSSLSETFDDGTQTLQGELTLALDKLAKNPSNPQLLAHEYQSK LSEYTLYRNAQSNTVKVIKDVDAIIQNFR
Research Area	Others
Source	E.coli
Target Names	mxiH
Protein Names	Recommended name: Protein mxiH
Expression Region	1-83aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	16.7 kDa
Protein Length	Full Length

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.