




Recombinant *Saccharomyces cerevisiae* GTP-binding nuclear protein GSP1/CNR1 (GSP1)

Product Code	CSB-EP333740SVGb0
Abbreviation	Recombinant <i>Saccharomyces cerevisiae</i> GSP1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P32835
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	SAPAANGEVPTFKLVLVGDGGTGKTTFVKRHLTGEFEKKYIATIGVEVHPLSFY TNFGEIKFDVWDTAGQEKFGLRDGYINAQCAIIMFDVTSRITYKNVPNWHR DLVRVCENIPIVLCGNKVDVKERKVKAKTITFHRKKNLQYYDISAKSNYNFEKPF LWLARKLAGNPQLEFVASPALAPPEVQVDEQLMQQYQQEMEQATALPLPDED DADL
Research Area	Cell Biology
Source	E.coli
Target Names	GSP1
Expression Region	2-219aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged
Mol. Weight	30.7 kDa
Protein Length	Full Length of Mature Protein

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.