




Recombinant Human O-acetyl-ADP-ribose deacetylase 1 (OARD1)

Product Code	CSB-EP896727HU
Abbreviation	Recombinant Human OARD1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9Y530
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	ASSLNEDPEGSRITYVKGDLFACPKTDSLHCISEDCRMGAGIAVLFFKKKFGGV QELLNQKKSGEVAVLKRDGRYIYYLITKKRASHKPTYENLQKSLEAMKSHCL KNGVTDLSMPRIGCGLDRLQWENVSAMIEEVFEATDIKITVYTL
Research Area	Others
Source	E.coli
Target Names	OARD1
Protein Names	Recommended name: O-acetyl-ADP-ribose deacetylase C6orf130 EC= 3.5.1.-
Expression Region	2-152aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	24.3 kDa
Protein Length	Full Length of Mature Protein

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.