




Recombinant Mouse Interleukin-1 receptor-associated kinase 3 (Irak3)

Product Code	CSB-EP011811MO
Abbreviation	Recombinant Mouse Irak3 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8K4B2
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Mus musculus (Mouse)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	MAGRCGARGALSPQLLLFDLPPALLGELCGILDSCDGPLGWRGLAERLSNSW LDVRHIEKYVNQGKSGTRELLWSWAQKNKTIGDLLEVLQDMGHQRAIHLIINYG VSWTPSVQTHHELPFSPFPPEVKHACRENDPGPLEPANVTVDNVLVPEHNEK GTLQKTPISFQSILEGTKHFHKDFLIGEGEIFEVYRVDIRNQAYAVKLFKQEKKM QLKKHWKRFLSELEVLLLFRHPHILELAAYFTETEKLCLVYPYMSNGTLFDRLQ CTNGTTPLSWHVRISVLIGIAKAIQYLHNTQPCAVICGNVSSANILLDDQLQPKL TDFAAAHFRPNLEQQSSTINMTGGGRKHLWYMPEEYIRQGRLSVKTDVYSFGI VIMEVLTGCKVVLDDPKHVQLRDLLMELMEKRGDLSCLFDRKIPPCPRNFS AKLFSLAGRCVATKAKLRPTMDEVLSSLESTQPSLYFAEDPPTSLKSFRCPSPL FLDNVPSIPVEDDENQNNHSVPPKEVLGTDRTVQKTPFECSQSEVTFLGLDRN RGNRGSEADCNVPSSSHEECWSPELVAPSQDLSPTVISLGSSEWVPGHSHYG SKPMEKRCSSGLFCSEHEQSKKQ
Research Area	Immunology
Source	E.coli
Target Names	Irak3
Protein Names	Recommended name: Interleukin-1 receptor-associated kinase 3 Short name= IRAK-3 EC= 2.7.11.1 Alternative name(s): IL-1 receptor-associated kinase M Short name= IRAK-M
Expression Region	1-609aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	73.5 kDa
Protein Length	Full Length



Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.