




# Recombinant Rat 17-beta-hydroxysteroid dehydrogenase 13 (Hsd17b13)

<b>Product Code</b>	CSB-EP711723RA
<b>Abbreviation</b>	Recombinant Rat Hsd17b13 protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5M875
<b>Storage Buffer</b>	Tris-based buffer,50% glycerol
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Purity</b>	≥ 85% as determined by SDS-PAGE.
<b>Sequence</b>	ESLVKFFIPQRRKSVAGQTVLITGAGHGIGRLTAYEFAKQKSRLVLWDISKHGV EETAACKRKLGAVVHVVFVDCSNRAEIYKSVDQVKKEVGDIEIVVNNAGAIYPA DLLSTKDEEITKTFEVNILGHFWIIKALLPSMLRRNSGHIVTVASVCGHRVIPYLIP YCSSKFAAVGFHRALTAELDTLGKTGIKTSCLCPVVFVNTGFTKNPSTRLWPVLE PDEVARSLIDGILTNNKMIFVPSYINISLIVEMFFPERVLKAINRIQNIQFEAIVGHR TKRK
<b>Research Area</b>	Others
<b>Source</b>	E.coli
<b>Target Names</b>	Hsd17b13
<b>Protein Names</b>	Recommended name: 17-beta-hydroxysteroid dehydrogenase 13 Short name= 17-beta-HSD 13 EC= 1.1.-.- Alternative name(s): Short-chain dehydrogenase/reductase 9
<b>Expression Region</b>	20-300aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 10xHis-tagged
<b>Mol. Weight</b>	37.4 kDa
<b>Protein Length</b>	Full Length of Mature Protein
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.