



Recombinant Human Tripartite motif-containing protein 72 (TRIM72), partial

Product Code	CSB-EP744394HU1(M)
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6ZMU5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	SAAPGLLHQELSCPLCLQLFDAPVTAECGHSFCRACLGRVAGEPAADGTVLC PCCQAPTRPQALSTNLQLARLVEGLAQVPQGHCEEHLDP LSIYCEQDRALVC GVCASLGSHRGHRLLPAAEAHARLKTQLPQQKLQLQEACMRKEKSVAVLEHQ LVEVEETVRQFRGAVGEQLGKMRVFLAALEGLDREAERVGEAGVALRREL GSLNSYLEQLRQMEKVL EEVADKPQTEFLMKYCLVTSRLQKILAEAPPPARLDI QLPIISDDFKFQVWRKMFALMPALEELTFDPSSAHP SLVVSSSGRRVECSEQ KAPPAGEDPRQFDKAVAVVAHQQLSEGEHYWEVDVGD KPRWALGVIAAEAP RRGRLHAVPSQGLWLLGLREGKILEAHVEAKEPRALRSPERRPTRIGLYLSFG DGVLSFYDASDADALVPLFAFHERLPRPVY PFFDVCWHDKGKNAQPLLLVGP EGAEA
Research Area	Cell Biology
Source	E.coli
Target Names	TRIM72
Expression Region	2-477aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal GST-tagged
Mol. Weight	79.4 kDa
Protein Length	Partial
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.