




Recombinant Pig Interleukin-4 receptor subunit alpha (IL4R), partial

Product Code	CSB-YP771258PI
Abbreviation	Recombinant Pig IL4R protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q863Z5
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Sus scrofa (Pig)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	VRVLEWPICLSDYVSTSTCEWRMAGPVNCSAEFRLSYQLKFFNTENHTTCVP ENRAGSVCVCHMLMESIVIVDTYQLDLWAGEQLLWNSSFKPSQNVKPLAPRN LMVHANISHTWLLTWSNPYPSESYLYSELTLYLVNISNENDPTDFRIYNVTYLG TLRFPANTLKSGAAYSARVKAWAQRYNSTWSEWSPSVKWLNYEOPLEQR
Research Area	Cancer
Source	Yeast
Target Names	IL4R
Protein Names	Recommended name: Interleukin-4 receptor subunit alpha Short name= IL-4 receptor subunit alpha Short name= IL-4R subunit alpha Short name= IL-4R-alpha Short name= IL-4RA Alternative name(s): CD_antigen= CD124
Expression Region	33-240aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	26.2
Protein Length	Extracellular Domain
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.