




Recombinant Human Glucagon receptor (GCGR), partial

Product Code	CSB-EP009316HU1
Abbreviation	Recombinant Human GCGR protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P47871
Storage Buffer	Tris-based buffer, 50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	AQVMDFLFEKWKLYGDQCHHNLSELLPPPTLVNRTFDKYSCWPDPNTTANISCPWYLPWVHHKVQHRFVFKRCGPDGQWVRGPRGQPWRDASQCQMDGEIEVQKEVAK
Research Area	Cancer
Source	E.coli
Target Names	GCGR
Protein Names	Recommended name: Glucagon receptor Short name= GL-R
Expression Region	26-136aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	20.5 kDa
Protein Length	Partial

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.