




Recombinant Human Gastric inhibitory polypeptide receptor (GIPR), partial

Product Code	CSB-EP009438HU1
Abbreviation	Recombinant Human GIPR protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P48546
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	RAETGSKGQTAGELYQRWERYRRECQETLAAAEPSPGLACNGSFDMYVCW DYAAPNATARASCPWYLPWHHHVAAGFVLRQCGSDGQWGLWRDHTQCENP EKNEAFLDQRLILERLQ
Research Area	Others
Source	E.coli
Target Names	GIPR
Protein Names	Recommended name: Gastric inhibitory polypeptide receptor Short name= GIP-R Alternative name(s): Glucose-dependent insulinotropic polypeptide receptor
Expression Region	22-138aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	17.5 kDa
Protein Length	Partial

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.