




Recombinant Rotavirus A Non-structural glycoprotein 4, partial

Product Code	CSB-EP889507RFU
Abbreviation	Recombinant Rotavirus A Non-structural glycoprotein 4 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9PYC8
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Rotavirus A (isolate RVA/Cow/United States/B223/1983/G10P8[11]) (RV-A)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	PTMKIALKTSKCSYKVVKYCIVTIFNTLLKLAGYKEQITTKDEIEKQMERVVKEM RRHFKMIDKLTTREIEQVGLLKRIHDKLDIRAVDEIDMTKEINQKNVRTLEEW EWGKNPYEPKEVTAAM
Research Area	Others
Source	E.coli
Target Names	N/A
Protein Names	Recommended name: Non-structural glycoprotein 4 Short name= NSP4 Alternative name(s): NCVP5 NS28
Expression Region	52-175aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-KSI-tagged
Mol. Weight	30.1 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.