




Recombinant Middle East respiratory syndrome-related coronavirus Orf1ab, partial

| | |
|--------------------------|--|
| Product Code | CSB-EP2883GJE |
| Abbreviation | Recombinant Middle East respiratory syndrome-related coronavirus Orf1ab protein, partial |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | A0A0G3U5R1 |
| Storage Buffer | Tris-based buffer,50% glycerol |
| Product Type | Recombinant Proteins |
| Immunogen Species | Middle East respiratory syndrome-related coronavirus (MERS-CoV) |
| Purity | Greater than 90% as determined by SDS-PAGE. |
| Sequence | SKDSNFLNRVRSIVNARIEPCSSGLSTDVVFRAFDICNYKAKVAGIGKYYKTN TCRFVELDDQGHHLDSYFVVKRHTMENYELEKHCYDLLRDCDAVAPHDFFIFD VDKVKTPHIVRQRLTEYTMMDLVYALRHFDQNSEVLKAILVKYGCCDVTYFEN KLWFDFVENPSVIGVYHKLGERVRQAILNTVKFCDHMKAGLVGVLTLDNQDL NGKWYDFGDFVITQPGSGVAIVDSYYSYLMPVLSMTDCLAAETHRDCDFNKP LIEWPLTEYDFTDYKVQLFEKYFKYWDQTYHANCVNCTDDRCVLHCANFNVL AMTMPKTCFGPIVRKIFVDGVPFVSCGYHYKELGLVMNMDVSLHRHRLSLKE LMMYAADPAMHIASSNAFLDLRTSCFSVAALTTGLTFQTVRPGNFNQDFYDFV VSKGFFKEGSSVTLKHFFFAQDGNAAITDYNYSYNLPTMCDIKQMLFCMEVV NKYFEIYDGGCLNASEVVNNLDKSAGHPFNKFGKARVYYESMSYQEQLDF AMTKRNVIPMTQMNLYAISAKNRARTVAGVSILSTMTNRQYHQMLKSMAA TRGATCVIGTTKFGYGGWDFMLKTLTKYKDVNDPHLMGWDYPKCDRAMPNMCRI FASLILARKHGTCCTTRDRFYRLANECAQVLSEYVLCGGGGYVVKPGGTSSGD ATTAYANSVFNILQATTANVSALMGANGNKIVDKEVKDMQFDLYVNVYRSTSP DPKFVDKYYAFLNKHFSMMILSDDGVVCYNSDYAAKGYIAGIQNFKETLYYQN NVFMSEAKCWVETDLKKGPFHEFCSQHTLYIKDGGDGYFLPYDPSPRILSAGCF VDDIVKTDGTLMVERFVSLAIDAYPLTKHEDIEYQNVFWVYLQYIEKLYKDLTGH MLDSYSVMLCGDNSAKFWEEAFYRDLYSSPTTLQ |
| Research Area | Others |
| Source | E.coli |
| Expression Region | 4378-5310aa |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag-Free |
| Mol. Weight | 106.9 kDa |



Protein Length

Partial

Image

(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.