




Recombinant Human Placenta-expressed transcript 1 protein (PLET1)

Product Code	CSB-EP754399HU
Abbreviation	Recombinant Human PLET1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6UQ28
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	TFIRYSSTCFTFDEYYTITLDIKASSHIYESNAVYSVFVPVND SVYAVVMKTLDE NSDSAGLWQRADKNCYSNSTYYVKDQYMTVLEAQWQAPEPENITEVEIQAFT VQIRALPILSTLKLREKLSTLALAAKIPQSSAFKPFFMITPKSIRLEGLANQVFS
Research Area	Cell Biology
Source	E.coli
Target Names	PLET1
Protein Names	Recommended name: Placenta-expressed transcript 1 protein
Expression Region	26-187aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	24.4 kDa
Protein Length	Full length of the mature protein

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.