




# Recombinant Human Nuclear receptor ROR-alpha (RORA), Partial

<b>Product Code</b>	CSB-YP020069HU1(F1)
<b>Abbreviation</b>	Recombinant Human RORA protein, partial
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P35398
<b>Storage Buffer</b>	Tris-based buffer, 50% glycerol
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	Greater than 85% as determined by SDS-PAGE.
<b>Sequence</b>	GETSPTVSMAELEHLAQNISKSHLETQCYLREELQQITWQTFLQEEIENYQNK QREVMWQLCAIKITEAIQYVVEFAKRIDGFMELCQNDQIVLLKAGSLEVVFIRM CRAFDSQNNTVYFDGKYASPDVFKSLGCEDFISVFVFEFGKSLCSMHLTEDEIAL FSAFVLMASADRSWLQEKVKIEKLQKIQALQHVLPQKNHREDGILTKLICKVSTL RALCGRHTEKLMFAFKAIYPDIVRLHFPPLYKELFTSEFEPAMQIDG
<b>Research Area</b>	Cancer
<b>Source</b>	Yeast
<b>Target Names</b>	RORA
<b>Protein Names</b>	Recommended name: Nuclear receptor ROR-alpha Alternative name(s): Nuclear receptor RZR-alpha Nuclear receptor subfamily 1 group F member 1 Retinoid-related orphan receptor-alpha
<b>Expression Region</b>	262-523aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-tagged
<b>Mol. Weight</b>	35.4
<b>Protein Length</b>	Partial
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.