




Recombinant Human Melanoma-associated antigen 1 (MAGEA1)

Product Code	CSB-MP013323HUd9
Abbreviation	Recombinant Human MAGEA1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P43355
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MSLEQRSLHCKPEEALAEQAQAEALGLVCVQAATSSSSPLVLGTLEEVPTAGSTD PPQSPQGASAFPTTINFTRQRQPSEGSSSREEEGPSTSCILESLFRAVITKKVA DLVGFLLLLKYRAREPVTKAEMLESVIKNYKHCFPEIFGKASESLQLVFGIDVKEA DPTGHSYVLVTCLGLSYDGLLDGNQIMPKTGFLIIVLMIAMEGGHAPEEEIWE ELSVMEVYDGREHSAYGEPKLLTQDLVQEKYLEYRQVPDSDPARYEFLWGP RALAETSYVKVLEYVIKVSARVRFFFPSLREAAALREEEEGV
Research Area	Cancer
Source	Mammalian cell
Target Names	MAGEA1
Expression Region	1-309aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	C-terminal hFc1-tagged
Mol. Weight	63.2 kDa
Protein Length	Full Length

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.