




Recombinant Mouse C-type lectin domain family 4 member K (Cd207), partial

Product Code	CSB-MP837695MO
Abbreviation	Recombinant Mouse Cd207 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8VBX4
Product Type	Recombinant Proteins
Immunogen Species	Mus musculus (Mouse)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	QAVFYPRLMGKILDVKSDAQMLKGRVDNISTLGSDLKTERGRVDDAEVQMQUIV NTTLKRVRSQLSLETSMKIANDQLQILMSWGEVDSLAKIPELKRDLDKASAL NTKVQGLQNSLENVNKLKQQSDILEMVARGWKYFSGNFYFSRTPKWTWYSA EQFCISRKAHLTSVSSESEKFLYKAADGIPHWIGLTAKAGSEGDWYWVDQTSF NKEQSRRFWIPGEPNAGNNEHCANIRVSALKCWNDGPCDNTFLFICKRPYV QTTE
Research Area	Metabolism
Source	Mammalian cell
Target Names	Cd207
Protein Names	Recommended name: C-type lectin domain family 4 member K Alternative name(s): Langerin CD_antigen= CD207
Expression Region	63-331aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged
Mol. Weight	34.8 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.