




Recombinant Mouse Complement C1q-like protein 3 (C1qI3)

Product Code	CSB-BP863674MOB1
Abbreviation	Recombinant Mouse C1qI3 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9ESN4
Product Type	Recombinant Proteins
Immunogen Species	Mus musculus (Mouse)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	HYEMLGTCRMVCDPYGGTKAPSTAATPDRGLMQSLPTFIQGPKEAGRPGK AGPRGPPGEPGPPGPVGPPEKGEPRQQLPGPPGAPGLNAAGAI SAATYS TVPKIAFYAGLKRQHEGYEV LKFDDVVTNLGNHYDPTTGKFTCSIPGIYFFTYH VLMRGGDGTSMWADLCKNNQVRASAI AQDADQNYDYASNSVVLHLEPGDEV YIKLDGGKAHGGNNNKYSTFSGFIIYAD
Research Area	Others
Source	Baculovirus
Target Names	C1qI3
Expression Region	21-255aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	28.6
Protein Length	Full Length of Mature Protein

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.