




# Recombinant Mouse Interleukin-1 receptor-associated kinase 3 (Irak3)

<b>Product Code</b>	CSB-BP011811MO
<b>Abbreviation</b>	Recombinant Mouse Irak3 protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8K4B2
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	Greater than 85% as determined by SDS-PAGE.
<b>Sequence</b>	MAGRCGARGALSPQLLLFDLPPALLGELCGILDSCDGPLGWRGLAERLSNSW LDVRHIEKYVNQGKSGTRELLWSWAQKNKTIGDLLEVLQDMGHQRAIHLIINYG VSWTPSVQTHHELPFSPFPPEVKHACRENDPGPLEPANVTVDNVLVPEHNEK GTLQKTPISFQSILEGTHKHFHKDFLIGEGEIFEVYRVDIRNQAYAVKLFKQEKKM QLKKHWKRFLSELEVLLLFRHPHILELAAYFTETEKLCLVYPYMSNGTLFDRLQ CTNGTTPLSWHVRISVLIGIAKAIQYLHNTQPCAVICGNVSSANILLDDQLQPKL TDFAAAHFRPNLEQQSSTINMTGGGRKHLWYMPEEYIRQGRLSVKTDVYSFGI VIMEVLTGCKVVLDDPKHVQLRDLLMELMEKRGDLSCLSFLDRKIPPCPRNFS AKLFSLAGRCVATKAKLRPTMDEVLSSESTQPSLYFAEDPPTSLKSFRCPSPL FLDNVPSIPVEDDENQNNHNSVPPKEVLGTDRTVQKTPFECSQSEVTFLGLDRN RGNRGSEADCNVPSSSHEECWSPELVAPSQDLSPTVISLGSSWEVPGHSYG SKPMEKRCSSGLFCSEHEQSKKQ
<b>Research Area</b>	Immunology
<b>Source</b>	Baculovirus
<b>Target Names</b>	Irak3
<b>Protein Names</b>	Recommended name: Interleukin-1 receptor-associated kinase 3 Short name= IRAK-3 EC= 2.7.11.1 Alternative name(s): IL-1 receptor-associated kinase M Short name= IRAK-M
<b>Expression Region</b>	1-609aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 10xHis-tagged and C-terminal Myc-tagged
<b>Mol. Weight</b>	72.3
<b>Protein Length</b>	Full Length
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.