




# Recombinant Human Mothers against decapentaplegic homolog 6 (SMAD6)

<b>Product Code</b>	CSB-BP021792HU
<b>Abbreviation</b>	Recombinant Human SMAD6 protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	O43541
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	Greater than 85% as determined by SDS-PAGE.
<b>Sequence</b>	MFRSKRSGLVRRLLWRSRVVPDREEGGSGGGGGDEDEDGSLGSRAEPAPRAR EGGGCGRSEVRPVAPRRPRDAVGQRGAQGAGRRRRAGGPPRPMSEPGAG AGSSLLDVAEPGGPGWLPESDCETVTCCLFSERDAAGAPRDASDPLAGAALE PAGGGRSREARSLLLLLEQELKTVTYSLKRLKERSLDTLLEAVESRGGVPGG CVLVPRADLRLGGQPAPPQLLLGRLFRWPDQLHAVELKPLCGCHSFAAAADG PTVCCNPYHFSRLCGPESPPPPYSRSLSPRDEYKPLDLS DSTLSYTETETATNSLI TAPGEFSDASMSPDATKPSHWCSVAYWEHRTRVGRLYAVYDQAVSIFYDLPQ GSGFCLGQLNLEQRSESVRRTRSKIGFGILLSKEPDGVWAYNRGEHPIFVNSP TLDAPGGRALVVRKVPPGYSIKVFDFERSGLQHAPEPDAADGPYDPNSVRISF AKGWGPCYSRQFITSCPCWLEILLNNPR
<b>Research Area</b>	Epigenetics and Nuclear Signaling
<b>Source</b>	Baculovirus
<b>Target Names</b>	SMAD6
<b>Protein Names</b>	Recommended name: Mothers against decapentaplegic homolog 6 Short name= MAD homolog 6 Short name= Mothers against DPP homolog 6 Alternative name(s): SMAD family member 6 Short name= SMAD 6 Short name= Smad6 Short na
<b>Expression Region</b>	1-496aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 10xHis-tagged and C-terminal Myc-tagged
<b>Mol. Weight</b>	57.4
<b>Protein Length</b>	Full Length
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.