




Recombinant Mouse Desmoglein-3 (Dsg3), partial

Product Code	CSB-BP007205MO
Abbreviation	Recombinant Mouse Dsg3 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O35902
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	EWVKFAKPCREREDNSRRNPIAKITSDFQKNQKITYRISGVGIDQPPFGIFVVD PNNGDINITAIVDREETPSFLITCRALNALGQDVERPLILTVKILDVNDNPPIFSQT IFKGEIEENSASNSLVMILNATDADEPNHMNSKIAFKIVSQEPAGMSMFLISRNT GEVRTLTSSLDREQISSYHLVSGADNDGTGLSTQCECSIKIKDVNDNFPVLR SQYSARIEENTLNAELLRFQVTDWDEEYTDNWLAVYFFTSGNEGNWFEIETDP RTNEGILKVVKALDYEQVQSMQFSIAVRNKAEFHQSVISQYRVQSTPVTIQVID VREGISFRPPSKTFTVQRGVSTNKLVGILGTYQATDEDTGKAASSVRYVLGR NDGGLLVIDSKTAQIKFVKNIDRDSTFIVNKTISAEVLAIDENTGKTSTGTIYVEVP SFNENCPSVVLEKKDICTSSPSVTLVRLDRGKYTGPYTVSLEEQLKLPVM WTITTLNATSALLQAQQQVSPGVYNVPVIVKDNQDGLCDTPESLTLTVCQCDD RSMCRAPIPSREPNTYGESSWRLGP
Research Area	Others
Source	Baculovirus
Target Names	Dsg3
Protein Names	Recommended name: Desmoglein-3 Alternative name(s): 130 kDa pemphigus vulgaris antigen homolog
Expression Region	50-617aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	64.3
Protein Length	Partial

Image



Based on the SEQUEST from database of Baculovirus host and target protein, the LC-MS/MS Analysis result of CSB-BP007205MO could indicate that this peptide derived from Baculovirus-expressed Mus musculus (Mouse) Dsg3.

