




Recombinant Vaccinia virus Protein OPG161 (OPG161), partial

Product Code	CSB-BP300755VAA1
Abbreviation	Recombinant Vaccinia virus OPG161 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P68616
Product Type	Recombinant Proteins
Immunogen Species	Vaccinia virus (strain Copenhagen) (VACV)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	VRLNQCMSANEAAITDAAVAVAAASSTHRKVASSTTQYDHKESCNGLYYQGS CYILHSDYQLFSDAKANCTAESSTLPNKSDVLITWLIDYVEDTWGSDGNPITKT TSDYQDSDVSQEVRKYFCVKTMN
Research Area	Others
Source	Baculovirus
Target Names	OPG161
Protein Names	Recommended name: Protein A33
Expression Region	57-185aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	18.3
Protein Length	Partial

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.