



SAT1 Antibody, Biotin conjugated

Product Code	CSB-PA02759D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21673
Immunogen	Recombinant Human Diamine acetyltransferase 1 protein (5-171AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Diamine acetyltransferase 1 (EC 2.3.1.57) (Polyamine N-acetyltransferase 1) (Putrescine acetyltransferase) (Spermidine/spermine N(1)-acetyltransferase 1) (SSAT) (SSAT-1), SAT1, SAT
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SAT1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.