



Defensin D2 Antibody, HRP conjugated

Product Code	CSB-PA305599LB01FKI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P81571
Immunogen	Recombinant Spinacia oleracea Defensin D2 protein (1-52AA)
Raised In	Rabbit
Species Reactivity	Spinacia oleracea
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Defensin D2 (Antimicrobial peptide D2) (So-D2) (Fragment)
Immunogen Species	Spinacia oleracea
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.