




Recombinant Gorilla gorilla gorilla Fibroblast growth factor

Product Code	CSB-CF2954GKM
Abbreviation	Recombinant Gorilla gorilla gorilla Fibroblast growth factor protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	G3QFY8
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose
Product Type	Transmembrane Proteins
Immunogen Species	Gorilla gorilla gorilla (Western lowland gorilla)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	HPIPDSSPLLQFGGQVRQRYLYTDDAQQTEAHLEIREDGTVGGAADQSPESLL QLKALKPGVIQILGVKTSRFLCQRPDGLYGSLSHFDPEACSFRELLLEDGYNVY QSEAHGLPLHLPGNKSPHRDPAPRGPARFLPLPGLPPAPPEPPGILAPQPPDV GSSDPLSMVGPSQGRSPSYAS
Research Area	Others
Source	in vitro E.coli expression system
Target Names	FGF21
Protein Names	Recombinant Gorilla gorilla gorilla Fibroblast growth factor(FGF21)
Expression Region	29-209aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged
Mol. Weight	22.2 kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.