




Recombinant Human Transforming growth factor beta-2 (TGFB2), partial

Product Code	CSB-EP023448HU1
Relevance	TGF-beta 2 has suppressive effects on interleukin-2 dependent T-cell growth.
Abbreviation	Recombinant Human TGFB2 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P61812
Alias	BSC-1 cell growth inhibitorCetermin;Glioblastoma-derived T-cell suppressor factor ;G-TSFPolyergin
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	LDAAYCFRNVQDNCCLRPLYIDFKRDLGWKWIHEPKGYNANFCAGACPYLWS SDTQHRSRVLSLYNTINPEASASPCCVSQDLEPLTILYYIGKTPKIEQLSNMIVKS CKC
Source	E.coli
Target Names	TGFB2
Expression Region	304-413aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	16.7kDa
Protein Length	Partial

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

