





Recombinant Human Carbonic anhydrase 9 (CA9), partial

Product Code	CSB-EP614990HU
Abbreviation	Recombinant Human CA9 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q16790
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Source	E.coli
Target Names	CA9
Protein Names	Recommended name: Carbonic anhydrase 9 EC= 4.2.1.1 Alternative name(s): Carbonate dehydratase IX Carbonic anhydrase IX Short name= CA-IX Short name= CAIX Membrane antigen MN P54/58N Renal cell carcinoma-associated antigen
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal GST-tagged
Protein Length	Partial


Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Based on the SEQUEST from database of E.coli host and target protein, the LC-MS/MS Analysis result of CSB-EP614990HU could indicate that this peptide derived from E.coli-expressed Homo sapiens (Human) CA9.



Based on the SEQUEST from database of E.coli host and target protein, the LC-MS/MS Analysis result of CSB-EP614990HU could indicate that this peptide derived from E.coli-expressed Homo sapiens (Human) CA9.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.