




Recombinant Human T-cell leukemia virus 1 HTLV-1 basic zipper factor (HBZ)

Product Code	CSB-EP313738HQJ
Abbreviation	Recombinant Human T-cell leukemia virus 1 HBZ protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0C746
Product Type	Recombinant Proteins
Immunogen Species	Human T-cell leukemia virus 1 (strain Japan ATK-1 subtype A) (HTLV-1)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MVNFVSAGLFRCLPVSCPEDLLVEELVDGLLSLEEELKDKEEEEAVLDGLLSLEEESRGRLRRGPPGEKAPPRGETHRDRQRRAEKRRKRKKEREKEEEKQTAEYLKRKEEEKARRRRRAEKKAADVARRKQEEQERRERKWRQGAEKAKQHSARK EKMQELGIDGYTRQLEGEVESLEAERRKLLQEKEDLMGEVNYWQGRLEAMWLQ
Research Area	Cardiovascular
Source	E.coli
Target Names	HBZ
Protein Names	Recommended name: HTLV-1 basic zipper factor Short name= HBZ
Expression Region	1-209aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	32.3 kDa
Protein Length	Full length

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.