




Recombinant Human Atrial natriuretic peptide receptor 3 (NPR3), partial

Product Code	CSB-EP016025HU1
Abbreviation	Recombinant Human NPR3 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P17342
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	GGVGGGGGGGAGIGGGRQEREALPPQKIEVLVLLPQDDSYLFSLTRVRPAIEYA LRSVEGNGTGRLLPPGTRFQVAYEDSDCGNRALFSLVDRVAAARGAKPDLIL GPVCEYAAAPVARLASHWDLPLMSAGALAAGFQHKDSEYSHLTRVAPAYAKM GEMMLALFRHHHWSRAALVYSDDKLERNICYFTLEGVHEVFQEEGLHTSIYSF DETKDLLEDIVRNIQASERVVIMCASSDTIRSIMLVHRHGMTSGDYAFFNIEL FNSSSYGDGSWKRGDKHDFEAKQAYSSLQTVTLRTVKPEFEKFSMEVKSSV EKQGLNMEDYVNMFVEGFHDAILLYVLALHEVLRAGYSKKDGGKIIQQTWNRT FEGIAGQVSIANGDRYGD FSVIAMTDVEAGTQEVIQDYFGKEGRFEMRPNVK YPWGPLKLRIDENRIVEHTNSSPCKSSGGLEE
Research Area	Signal Transduction
Source	E.coli
Target Names	NPR3
Protein Names	Recommended name: Atrial natriuretic peptide receptor 3 Alternative name(s): Atrial natriuretic peptide clearance receptor Atrial natriuretic peptide receptor type C Short name= ANP-C Short name= ANPR-C Short name= NPR-C
Expression Region	27-481aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged
Mol. Weight	54.0 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.