




# Recombinant Human Hyaluronan and proteoglycan link protein 2 (HAPLN2)

<b>Product Code</b>	CSB-EP010131HU
<b>Abbreviation</b>	Recombinant Human HAPLN2 protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q9GZV7
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	Greater than 85% as determined by SDS-PAGE.
<b>Sequence</b>	DPASHPGPHYLLPPIHEVIHSHRGATATLPCVLGTTPPSYKVRWSKVEPGELR ETLILITNGLHARGYGPLGGRARMRRGHRLDASLVIAGVRLEDEGRYRCELING IEDESVALTSLLEGVVFPYQPSRGRYQFNYYEAKQACEEQDGRLATYSQLYQA WTEGLDWCNAGWLLEGSVRYPVLTARAPCGGRGRPGIRSYGPRDRMRDRY DAFCFTSALAGQVFFVPGRLTLSEAHAAACRRRGAVVAKVGHLYAAWKFSGLD QCDGGWLADGSVRFPIITTPRPRCGGLPDPGVRSFGFPRPQQAAYGTYCYAE N
<b>Research Area</b>	Neuroscience
<b>Source</b>	E.coli
<b>Target Names</b>	HAPLN2
<b>Protein Names</b>	Recommended name: Hyaluronan and proteoglycan link protein 2 Alternative name(s): Brain link protein 1
<b>Expression Region</b>	27-340aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 10xHis-tagged
<b>Mol. Weight</b>	40.8 kDa
<b>Protein Length</b>	Full Length of Mature Protein
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.