




Recombinant Mycobacterium smegmatis Porin MspA (mspA)

Product Code	CSB-EP369785MVX
Abbreviation	Recombinant Mycobacterium smegmatis mspA protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	A0QR29
Product Type	Recombinant Proteins
Immunogen Species	Mycobacterium smegmatis (strain ATCC 700084 / mc(2)155)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	GLDNELSLVDGQDRRLTVQQWDTFLNGVFPLDRNRLTREWFHSGRAKYIVAG PGADEFEGTLELGYQIGFPWSLGVGINFSYTTPNILIDDGDITAPPFGLNSVITP NLFPGVSISADLGNPGIQEVATFSVDVSGAEGGVAVSNAHGTVTGAAGGVLL RPFARLIASGTGDSVTTYGEPWNMN
Research Area	Others
Source	E.coli
Target Names	mspA
Protein Names	Recommended name: Porin MspA
Expression Region	28-211aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged
Mol. Weight	25.4 kDa
Protein Length	Full Length of Mature Protein

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.