




# Recombinant Rotavirus A Intermediate capsid protein VP6

<b>Product Code</b>	CSB-EP365972RIX
<b>Abbreviation</b>	Recombinant Rotavirus A Intermediate capsid protein VP6 protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P03530
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Rotavirus A (strain RVA/Human/United States/Wa/1974/G1P1A[8]) (RV-A)
<b>Purity</b>	≥ 85% as determined by SDS-PAGE.
<b>Sequence</b>	MEVLYSLSKTLKDARDKIVEGTLYSNVSDLIQQFNQMIVTMNGNDFQTGGIGNL PVRNWTDFDGLLGTTLLNLDANYVENARTIIEYFIDFIDNVCMDEMARESQRNG VAPQSEALRKLKLAGIKFKRINFDNSSEYIENWNLQNRRTGFVVFHKPNIFPYS SFTLNRSQPMHDNLMGTMWLNAGSEIQVAGFDYSCAINAPANIQQFEHIVQLR RALTTATITLLPDAERFSFPRVINSADGATTWFFNPVILRPNNVEVEFLLNGQIIN TYQARFGTIIARNFDAIRLLFQLMRPPNMTPAVNALFPQAQPFQHHATVGLTLR IESAVCESVLADANETLLANVTAVRQEYAI PVGPFVPPGMNWTEELITNYSRE DNLQRVFTVASIRSMLIK
<b>Research Area</b>	Others
<b>Source</b>	E.coli
<b>Target Names</b>	N/A
<b>Protein Names</b>	Recommended name: Intermediate capsid protein VP6
<b>Expression Region</b>	1-397aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-tagged
<b>Mol. Weight</b>	50.9 kDa
<b>Protein Length</b>	Full Length
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.