




Recombinant Bovine Ubiquitin-like protein ISG15 (ISG15)

Product Code	CSB-EP011843BO
Abbreviation	Recombinant Bovine ISG15 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O02741
Product Type	Recombinant Proteins
Immunogen Species	Bos taurus (Bovine)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	MGGDLTVKMLGGQEILVPLRDSMTVSELKQFIAQKINVPFAFQQRLAHLDSREV LQEGVPLVLQGLRAGSTVLLVQNCISILVRNDKGRSSPYEVQLKQTVVELKQ QVCQKERVQADQFWLSFEGRPMDDHEHPLEEYGLMKGCTVFMNLRRLRGG
Research Area	Cell Biology
Source	E.coli
Target Names	ISG15
Protein Names	Recommended name: Ubiquitin-like protein ISG15 Alternative name(s): Interferon-stimulated gene product 17 Ubiquitin cross-reactive protein Short name= BoUCRP
Expression Region	1-154aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	21.4 kDa
Protein Length	Full Length

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.