




Recombinant Nostoc sp. All1616 protein (all1616)

Product Code	CSB-EP2964FUZa0
Abbreviation	Recombinant Nostoc sp. All1616 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8YWJ7
Product Type	Recombinant Proteins
Immunogen Species	Nostoc sp. (strain PCC 7120 / SAG 25.82 / UTEX 2576)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	MRWVDDSYLLPAAKGLMSTSVNTYAKLAKALLINYSFDLSGYHVDELVNRW QKQYPADWLHLAVIEALYQGRYKAISVQQLLAFWQRRGQEIYHFNMEFERLIC SKFPESLTPMAASEQYSRQGKNQNQTLQLMSFKQQEQVKEEEEEPPTEKMLAL SSTSVTASIEVSVSQQEDYLGQPFSLNPDISTKLLPISVTHPPIGQFTPQTSDRS ESFTSKLKAISNENS
Research Area	Others
Source	E.coli
Target Names	all1616
Expression Region	1-227aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	31.9 kDa
Protein Length	Full Length

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.