




Recombinant Human metapneumovirus

Nucleoprotein (N)

Product Code	CSB-EP751042HDAM
Abbreviation	Recombinant Human metapneumovirus N protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6WBA1
Product Type	Recombinant Proteins
Immunogen Species	Human metapneumovirus(strain CAN97-83)(HMPV)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MSLQGIHLSDLSYKHAILKESQYTIKRDVGTTTAVTPSSLQQEITLLCGEILYAKH ADYKYAAEIGIQYISTALGSERVQQILRNSGSEVQVVLTRTYSLGKVKNNKGED LQMLDIHGVEKSWVEEIDKEARKTMATLLKESSGNIPQNQRPSAPDTPILLCV GALIFTKLASTIEVGLETTVRRANRVLSDALKRYPRMDIPKIARSFYDLFEQKVY YRSLFIEYGKALGSSSTGSKAESLFVNIFMQAYGAGQTMLRWGVIARSSNNIM LGHVSVQAELKQVTEVYDLVREMGPESGLLHLRQSPKAGLLSLANCPNFASV VLGNASGLGIIGMYRGRVPNTELFSAAESYAKSLKESNKINFSSLGLTDEEKEA AEHFLNVSDDSQNDYE
Research Area	Signal Transduction
Source	E.coli
Target Names	N
Protein Names	Recommended name: Nucleoprotein Short name= Protein N Alternative name(s): Nucleocapsid protein
Expression Region	1-394aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	51.0 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.