




Recombinant Mouse F-box only protein 32 (Fbxo32)

Product Code	CSB-EP887410MO
Abbreviation	Recombinant Mouse Fbxo32 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9CPU7
Product Type	Recombinant Proteins
Immunogen Species	Mus musculus(Mouse)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MPFLGQDWRSPGQSWVKTADGWKRFLDEKSGSFVSDLSSYCNKEVYSKENL FSSLNYDVAKKRKKDIQNSKTKTQYFHQEKWIYVHKGSTKERHGYCTLGEAF NRLDFSTAILDSRRFNYYVRLLELIAKSQLTSLSGIAQKNFMNILEKVVLKVLED QQNIRLIRELLQTLYTSLCTLVQRVGKSVLVGNINMWWYRMETILHWQQQLNSI QISRPAFKGLTITDLPVCLQLNIMQRLSDGRDLVSLGQAAPDLHVLSEDRLWLK RLCQYHFSEKQIRKRLILSDKGQLDWKKMYFKLVRCYPRREQYGVTLQLCKH CHILSWKGTDHPCTANNPESCSVSLSPQDFINLFKF
Research Area	Cardiovascular
Source	E.coli
Target Names	Fbxo32
Protein Names	Recommended name: F-box only protein 32 Alternative name(s): Atrogin-1 Muscle atrophy F-box protein Short name= MAFbx
Expression Region	1-355aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	48.9 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.