




Recombinant Vespa mandarinia Vespakinin-M

Product Code	CSB-EP762592VBR
Abbreviation	Recombinant Vespa mandarinia Vespakinin-M protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7M3T3
Product Type	Recombinant Proteins
Immunogen Species	Vespa mandarinia(Asian giant hornet)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	GRPPGFSPFRID
Research Area	Others
Source	E.coli
Target Names	N/A
Protein Names	Recommended name: Vespakinin-M
Expression Region	1-12aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-KSI-tagged
Mol. Weight	16.7 kDa
Protein Length	Full Length

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.