




# Recombinant Trypanosoma cruzi Cruzipain

<b>Product Code</b>	CSB-EP333107TQY
<b>Abbreviation</b>	Recombinant Trypanosoma cruzi Cruzipain protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P25779
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Trypanosoma cruzi
<b>Purity</b>	Greater than 85% as determined by SDS-PAGE.
<b>Sequence</b>	APAAVDWRARGAVTAVKDGQCGSCWAFSAIGNVECQWFLAGHPLTNLSEQ MLVSCDKTDSGCSGGLMNNAFEWIVQENNGAVYTEDSYPYASGEGISPPCTT SGHTVGATITGHVELPQDEAQIAAWLAVNGPVAVAVDASSWMTYTGGVMTSC VSEQLDHGVLLVGYNDSAAVPYWIKNSTTQWGEEGYIRIAKGSNQCLVKEE ASSAVVGGPGPTPEPTTTTTTSAPGPSYFVQMSCTDAACIVGCENVTLPTG QCLLTSGVSAIVTCGAETLTTEEVFLTSTHCSGPSVRSSVPLNKCNRLLRGSVE FFCGSSSSGRLADVDRQRHQPYPYHSRHRRL
<b>Research Area</b>	Others
<b>Source</b>	E.coli
<b>Target Names</b>	N/A
<b>Protein Names</b>	Recommended name: Cruzipain EC= 3.4.22.51Alternative name(s): Cruzaine Major cysteine proteinase
<b>Expression Region</b>	123-467aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 10xHis-tagged and C-terminal Myc-tagged
<b>Mol. Weight</b>	43.9 kDa
<b>Protein Length</b>	Full Length of Mature Protein
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.