




# Recombinant Human Zinc finger and BTB domain-containing protein 32 (ZBTB32)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-EP897538HU                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Abbreviation</b>      | Recombinant Human ZBTB32 protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                                                                                                          |
| <b>Uniprot No.</b>       | Q9Y2Y4                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Product Type</b>      | Recombinant Proteins                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Immunogen Species</b> | Homo sapiens (Human)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Purity</b>            | ≥ 85% as determined by SDS-PAGE.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Sequence</b>          | MSLPPIRLPSPYGS DRLVQLAARLRPALCDTLITVGSQEFPAHSLVLAGVSQQ<br>L GRRGQWALGEGISPSTFAQLLN FVYGESVELQPGELRPLQEAAALGVQSLE<br>EACWRARGDRAKKPD PGLKKHQEEPEKPSRNP ER ELGDPGEKQKPEQVSRT<br>GGREQEMLHKHSPPRGRPEMAGATQEAQQEQTRSKEKRLQAPVGQRGADG<br>KHGVLTLWLRNPGGSEESLRKLPGPLPPAGSLQTSVTPRPSWAEAPWLVG<br>G QPALWSILLMPPRYGIPFYHSTPTTGAWQEVWREQRIPLSLNAPKGLWSQ<br>NQ LASSSPTPGSLPQGP AQLSPGEMEESDQGHTGALATCAGHEDKAGC<br>PPRPH PPPAPPARSRPYACSVCGKRFS LKHQMETHYRVHTGEKPFSCSLCP<br>QRSRDF SAMTKHLRTHGAAPYRCSLCGAGCPSLASMQA HMRGHS<br>PSQLPPGWTIRST FLYSSSRPSRPSTSPCCPSSSTT |
| <b>Research Area</b>     | Epigenetics and Nuclear Signaling                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Source</b>            | E.coli                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Target Names</b>      | ZBTB32                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Protein Names</b>     | Recommended name: Zinc finger and BTB domain-containing protein 32<br>Alternative name(s): FANCC-interacting protein Fanconi anemia zinc finger protein Testis zinc finger protein Zinc finger protein 538                                                                                                                                                                                                                                                                                                                                               |
| <b>Expression Region</b> | 1-487aa                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| <b>Tag Info</b>          | N-terminal 10xHis-tagged and C-terminal Myc-tagged                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Mol. Weight</b>       | 60.4 kDa                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Protein Length</b>    | Full Length                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Image</b>             |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.