




Recombinant Human metapneumovirus Phosphoprotein (P)

Product Code	CSB-EP804307HDAM
Abbreviation	Recombinant Human metapneumovirus Phosphoprotein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8B9Q8
Product Type	Recombinant Proteins
Immunogen Species	Human metapneumovirus(strain CAN97-83)(HMPV)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MSFPEGKDILFMGNEAAKLAEAFQKSLRKPSHKRSQSIIGEKVNTVSETLELPTI SRPTKPTILSEPKLAWTDKGGAIKTEAKQTIKVMDDPIEEEEFTEKRVLPSSDGKT PAEKKLKPSTNTKKKVSFTPNEPGKYTKLEKDALLSDNEEEDAESSILTFEE RDTSSLSIEARLESIEEKLSMILGLLRTLNIATAGPTAARDGIRDAMIGIREELIADI IKEAKGKAAEMMEEEMNQRTKIGNGSVKLTEKAKELNKIVEDESTSGESEEEEE ELKDTQENNQEDDIYQLIM
Research Area	Signal Transduction
Source	E.coli
Target Names	P
Protein Names	Recommended name: Phosphoprotein Short name= Protein P
Expression Region	1-294aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	40.2 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.