




Recombinant Mouse Pulmonary surfactant-associated protein C (Sftpc)

Product Code	CSB-EP021174MO
Abbreviation	Recombinant Mouse Sftpc protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P21841
Product Type	Recombinant Proteins
Immunogen Species	Mus musculus(Mouse)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	FRIPCCPVHLKRLIVVVVVVLVVVVIVGALLMGL
Research Area	Cancer
Source	E.coli
Target Names	Sftpc
Protein Names	Recommended name: Pulmonary surfactant-associated protein C Short name= SP-C Alternative name(s): Pulmonary surfactant-associated proteolipid SPL(Val) SP5
Expression Region	24-58aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-KSI-tagged
Mol. Weight	19.1 kDa
Protein Length	Full Length of Mature Protein

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.