




Recombinant Human Hyaluronan and proteoglycan link protein 3 (HAPLN3)

Product Code	CSB-EP010132HU
Abbreviation	Recombinant Human HAPLN3 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q96S86
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	LPFYNGFYYSNSANDQNLGNGHGHKDLLNGVKLVVETPEETLFTYQGASVILPC RYRYEPALVSPRRVRVKWWKLSENGAPEKDVLVAIGLRHRSFGDYQGRVHLR QDKEHDVSLEIQDLRLEDYGRYRCEVIDGLEDESGLVELELRGVVFPYQSPNG RYQFNFHEGQQVCAEQAAVVASFEQLFRAWEEGLDWCNAGWLQDATVQYPI MLPRQPCGGPGLAPGVRVSYGPRHRRLHRYDVFCFATALKGRVYYLEHPEKLT LTEAREACQEDDATIAKVGQLFAAWKFHGLDRCDAGWLADGVSRYPVVHPPH NCGPPEPGVRSFGFPDPQSRLYGVCYRQH
Research Area	Signal Transduction
Source	E.coli
Target Names	HAPLN3
Protein Names	Recommended name: Hyaluronan and proteoglycan link protein 3
Expression Region	18-360aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged
Mol. Weight	42.8 kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.